

**SUBMISSION CASE NARRATIVE
NDMA**

MAXXAM L.I.M.S. No. A203768

PROJECT: Applied P & Ch Laboratory NDMA Analysis

I. Receipt

Samples were received at Maxxam on February 8, 2002.
Samples were received in good condition.

II. Holding Times

- A. Sample preparation: all holding times were met.
- B. Sample analysis: all holding times were met.

III. Method

The method followed was Maxxam's in-house method for NDMA analysis, Entitled "EXTRACTION & ANALYSIS OF NITROSAMINES AND NDMA BY HRMS" SOP # TO.1021.04.

IV. Preparation

Sample preparation proceeded normally. Samples were extracted on February 11, 2002.

V. Analysis


Analysis proceeded normally. Samples were analyzed on February 14, 2002.

- A. Calibration: All criteria were met.
- B. Mass Resolution: All criteria met.

000001

- C. Method Blank: All acceptance criteria were met for method blank.
- D. Laboratory Control Spike: A LCS and LCSDUP were analyzed but found to be contaminated. The Method Blank and Samples were not affected, results were reported except for the LSC and LCSDUP as per client agreement.
- E. Matrix spike/Matrix spike duplicate: MS and MSD were analyzed with these samples and had a RPD of 12%.
- F. Surrogate Standards: All samples and QC samples met surrogate Standard criteria
- G. Samples: Sample analysis proceeded normally.
- H. Glass blank: All acceptance criteria for glass blank were met.

I certify that this data package is in compliance with the terms and conditions agreed to by the client and Maxxam Analytics Inc., both technically and for completeness, except for any conditions noted above. Release of the data contained in this hardcopy data package has been authorized by the HRMS Strategic Business Unit Operational Manager, as verified by the following signature.



AnnMarie Wright, B.Sc.
Laboratory Operations Manager

SUMMARY OF SAMPLES SUBMITTED-NDMA									
(YYYY/MM/DD)									
SOTA SAMPLE NO.	MAXXAM	DATE	DATE	DATE	DATE	DATE	DATE	DATE	ASSOCIATED
	L.I.M.S. ID	SAMPLED	RECEIVED	EXTRACTED	ANALYZED	QC LABEL			
ER-13	A203768-714490	2002/02/05	2002/02/08	2002/02/11	2002/02/14	20020211			
MW-4-2	A203768-714491	2002/02/05	2002/02/08	2002/02/11	2002/02/14	20020211			
MW-4-2 MS	A203768-714492	2002/02/05	2002/02/08	2002/02/11	2002/02/14	20020211			
MW-4-2 MSD	A203768-714493	2002/02/05	2002/02/08	2002/02/11	2002/02/14	20020211			

000003

Glossary of Definitions

NDMA	N-Nitrosodimethylamine
OPR	Ongoing Performance & Recovery Standard (Matrix spike)
PAR	Performance & Recovery Standard (Spiking Mixture)
IPR	Initial Performance & Recovery Standard (Matrix spike)
K-D	Kuderna-Danish concentrator; a device used to concentrate the analytes in a solvent
LIMS	Laboratory Information Management System
MISA	Municipal Industrial Strategy for Abatement
EPA	see USEPA
USEPA	United States Environmental Protection Agency
CEPA	Canadian Environmental Protection Agency
amp	ampere
cm	centimetre
g	gram
h	hour
ID	internal diameter
OD	outside diameter
In.	inch
L	litre
M	Molecular ion
min	minute
mL	millilitre
mm	millimetre
m/z	mass-to-charge ratio
N	Normal; gram molecular weight of solute divided by hydrogen equivalent of solute, per litre of solution
mg	milligram 10^{-3} g
μ g	microgram 10^{-6} g
ng	nanogram 10^{-9} g
pg	picogram 10^{-12} g
fg	femtogram 10^{-15} g
ppm	parts per million (mg/L, mg/kg)
ppb	parts per billion (μ g/L, μ g/kg)
ppt	parts per trillion (ng/L, ng/kg)
ppq	parts per quadrillion (pg/L, pg/kg)
v/v	volume per unit volume
w/v	weight per unit volume
DCM	Dichloromethane (Methylene Chloride)
PFK	Perfluorokerosene
HIRES	High Resolution
GC	Gas Chromatography

MS Mass Spectrometry
HRMS High Resolution Mass Spectrometry

Acceptance Criteria

Values used by the laboratory in order to determine that a process is in control.

Accuracy It is the degree of agreement of a measured value with the true or expected value of the quantity of concern.

Analyte A Nitrosodimethylamine and/or 1,4-Dioxane parameter tested by a method.

Blind Sample It is a sample submitted for analysis whose composition is known to the submitter but unknown to the analyst. A blind sample is used to test the proficiency of a measurement process.

Calibration Standard (CAL)

Consist of a set of solutions containing known amounts of native & carbon-13-labelled NDMA and/or 1,4-Dioxane. These solutions are used to establish the relationship between the parameter's concentration & MS detector response over the expected range of sample concentration.

Calibration Verification Material

Consists of a calibration standard solution of intermediate level concentration (e.g. CS3), used to assess whether the initial calibration is still valid.

Certified Reference Material

It is a stable, homogenous, and well characterized reference material, one or more of whose property values are certified by repetitive analysis by several operators & different methodologies in one or more qualified laboratories of known precision & accuracy. This material is used to assess the accuracy of a measurement process.

CAS# Chemical Abstracts Compound Registry Number.

Control Sample

It is a reference material of known composition that is analyzed concurrently with test samples to evaluate the accuracy and/or precision of a measurement process.

EDL Estimated detection limit or detection limit.

Glassware Proof Rinse

It is the composite final solvent rinse of each piece of glassware intended for use in processing a batch of samples. Proof rinse samples are analyzed before sample processing begins.

Instrument Detection Limit

It is the smallest concentration/amount of analyte, in a solution containing only the analyte(s) of interest, which produces an instrumental response that satisfies all analyte detection & identification criteria.

IS

Internal Standard, a deuterated or ^{13}C -labelled analyte that is added to a sample extract prior to instrument analysis.

Isomer

A member of a group of compounds that differ from each other only in terms of locations of a specified number of common substituent atoms, or groups of atoms, on the parent compound.

Method Blank Laboratory control sample using reagents, purified water, soil or relevant matrix known to be free of contaminants.

Method Detection Limit (MDL)

It is the smallest test sample concentration/amount of analyte that produces an instrumental response that satisfies all analyte detection & identification criteria when the sample is processed & analyzed according to the requirements of a specific test method. Reported MDL values reflect the composite effect of sample-related variables as well as method-related variables.

MSDS

Material Safety Data Sheet

NIOSH

National Institute of Occupational Safety & Health

Precision

It is the degree of agreement between the data generated from repetitive measurements under specified conditions. It is generally reported as the standard deviation (SD) or relative standard deviation (RSD).

%D

Percent Difference.

Quality Assurance (QA)

It is a system of activities whose purpose is to provide the producer or user of a product with the assurance that the product meets a defined standard of quality. The system consists of two separate but related activities, quality control & quality assessment.

Quality Control (QC)

It is the overall system of activities whose purpose is to control the quality of a product so that it meets the needs of users.

Recovery Standards

They are selected compounds that are added to sample extracts immediately before instrumental analysis so that surrogate (internal standard) recoveries can be calculated.

RPD (%) Relative Percent Difference.

Relative Retention Factor (RRF)

It is the quotient of a target analyte response factor (instrument response per unit weight) divided by the response factor (RF) for its corresponding labelled surrogate. An RRF value remains constant over the range of concentration for which instrument response is linear.

RSD Relative Standard Deviation.

SDS Soxhlet/Dean-Stark extractor, an extraction device applied to the extraction of solid & semi-solid materials.

Spiked blank Laboratory control sample that has been fortified with native analytes of interest.

Stock Solution A solution containing an analyte that is prepared using a reference material traceable to EPA, the National Institute of Science & Technology (NIST), or a source that will attest to the purity & authenticity of the reference material.

Surrogate A compound whose composition and chemical properties are nearly identical to those of target analytes, but which is distinguishable from target analytes by some means of detection (i.e. MS). These include deuterated or ¹³C-labelled analogues of the target analytes, which are added to the sample prior to extraction or clean-up steps.

Window Defining Mixture

It is a solution containing the earliest & latest eluting congeners within each homologous group of target analytes on a specified GC column.

SAMPLE DATA

000008

Lab Name Maxxam Analytics Inc.

Matrix (soil/water): water

Sample wt/vol: 990 (g/mL) mL

Level (low/med) low

% Moisture Not applicable Decanted (Y/N): N

Concentrated Extract Volume 1000 (uL)

Injection Volume 2 (uL)

Acid Wash Cleanup (Y/N): N pH Not analyzed

Lab Sample ID: A203768-714490

Project Name: JPL

Lab File ID: KR0752

Date Received: February 8, 2002

Date Extracted: February 11, 2002

Lab Batch: 20020211

Date Analyzed: February 14, 2002

Calib. Ref.: 20020214

Time Analyzed: 18:52

Dilution Factor: 1

CAS No.	Compound	Conc. (ug/L)	Qualifier	EDL (ug/L)	RL (ug/L)
62-75-9	NDMA	0.000800	J	0.000390	0.00200
	Surrogate	Recovery (%)	Acceptance Criteria (%)		
76523-40-5	D6-NDMA	16	10-85		

000009

MW-4-2

Lab Name Maxxam Analytics Inc.

Matrix (soil/water): water

Sample wt/vol: 970 (g/mL) mL

Level (low/med) low

% Moisture Not applicable Decanted (Y/N): N

Concentrated Extract Volume 1000 (uL)

Injection Volume 2 (uL)

Acid Wash Cleanup (Y/N): N pH Not analyzed

Lab Sample ID: A203768-714491

Project Name: JPL

Lab File ID: KR0753

Date Received: February 8, 2002

Date Extracted: February 11, 2002

Lab Batch: 20020211

Date Analyzed: February 14, 2002

Calib. Ref.: 20020214

Time Analyzed: 19:04

Dilution Factor: 1

CAS No.	Compound	Conc. (ug/L)	Qualifier	EDL (ug/L)	RL (ug/L)
62-75-9	NDMA	0.00148	J	0.000390	0.00200
	Surrogate	Recovery (%)	Acceptance Criteria (%)		
76523-40-5	D6-NDMA	15	10-85		

000011

Applied P & Ch Laboratory

13760 Magnolia Ave. Chino CA 91710

Tel: (909) 590-1828 Fax: (909) 590-1498

APCL Analytical Report

Submitted to:

SOTA Environmental

Attention: Yu Zeng

16835 W. Bernardo Dr, Ste. 212

San Diego CA 92127

Tel: (858)485-8100 Fax: (858)485-0812

Service ID #: 801-021400

Collected by:

Collected on: 01/29/02

Received: 01/29/02

Extracted: N/A

Tested: N/A

Reported: 02/27/02

Sample Description: Water

Project Description: 00HW019 JPL

Analysis of Water Samples

Component Analyzed	Method	Unit	PQL	Analysis Result		
				MW-13	MW-16	MW-16D
				02-01400-1	02-01400-2	02-01400-3

NITROSAMINES BY HRMS ^(a)

PQL: Practical Quantitation Limit. MDL: Method Detection Limit. CRDL: Contract Required Detection Limit

N.D.: Not Detected or less than the practical quantitation limit.


"-": Analysis is not required.

J: Reported between PQL and MDL.

Listed Dilution Factors (DF) are relative to the method default DF. All unlisted DFs are 1.0

^(a) Subcontracted to Maxxam Analytics Inc. See attached.

Respectfully submitted,



Dominic Lau

Laboratory Director

Applied P & Ch Laboratory

**SUBMISSION CASE NARRATIVE
NDMA**

MAXXAM L.I.M.S. No. A203024

PROJECT: Applied P & Ch Laboratory NDMA Analysis

I. Receipt

Samples were received at Maxxam on February 1, 2002.
Samples were received in good condition.

II. Holding Times

- A. Sample preparation: all holding times were met.
- B. Sample analysis: all holding times were met.

III. Method

The method followed was Maxxam's in-house method for NDMA analysis, Entitled "EXTRACTION & ANALYSIS OF NITROSAMINES AND NDMA BY HRMS" SOP # TO.1021.04.

IV. Preparation

Sample preparation proceeded normally. Samples were extracted on February 4, 2002.

V. Analysis


Analysis proceeded normally. Samples were analyzed on February 4, 2002.

- A. Calibration: All criteria were met.
- B. Mass Resolution: All criteria met.

000001

- C. Method Blank: All acceptance criteria were met for method blank.
- D. Laboratory Control Spike: A LCS and LCSDUP were analyzed and they had a RPD of 7%.
- E. Matrix spike/Matrix spike duplicate: MS and MSD were analyzed with these samples and they had a RPD of 1%.
- F. Surrogate Standards: All samples and QC samples met surrogate Standard criteria
- G. Samples: Sample analysis proceeded normally.
- H. Glass blank: All acceptance criteria for the glass blank were met.

I certify that this data package is in compliance with the terms and conditions agreed to by the client and Maxxam Analytics Inc., both technically and for completeness, except for any conditions noted above. Release of the data contained in this hardcopy data package has been authorized by the HRMS Facility Co-ordinator, as verified by the following signature.



AnnMarie Wright, B.Sc.
Laboratory Operations Manager

SUMMARY OF SAMPLES SUBMITTED-NDMA									
(YYYY/MM/DD)									
SOTA SAMPLE NO.	MAXXAM	DATE	DATE	DATE	DATE	DATE	DATE	DATE	ASSOCIATED
	L.I.M.S. ID	SAMPLED	RECEIVED	EXTRACTED	ANALYZED	QC LABEL			
MW-13	A203024-709615	2002/01/29	2002/02/01	2002/02/04	2002/02/04	20020204			
MW-13 MSD	A203024-709616	2002/01/29	2002/02/01	2002/02/04	2002/02/04	20020204			
MW-13 MS	A203024-709617	2002/01/29	2002/02/01	2002/02/04	2002/02/04	20020204			
MW-16	A203024-709618	2002/01/29	2002/02/01	2002/02/04	2002/02/04	20020204			
MW-16D	A203024-709619	2002/01/29	2002/02/01	2002/02/04	2002/02/04	20020204			

000003

Glossary of Definitions

NDMA	N-Nitrosodimethylamine
OPR	Ongoing Performance & Recovery Standard (Matrix spike)
PAR	Performance & Recovery Standard (Spiking Mixture)
IPR	Initial Performance & Recovery Standard (Matrix spike)
K-D	Kuderna-Danish concentrator; a device used to concentrate the analytes in a solvent
LIMS	Laboratory Information Management System
MISA	Municipal Industrial Strategy for Abatement
EPA	see USEPA
USEPA	United States Environmental Protection Agency
CEPA	Canadian Environmental Protection Agency
amp	ampere
cm	centimetre
g	gram
h	hour
ID	internal diameter
OD	outside diameter
In.	inch
L	litre
M	Molecular ion
min	minute
mL	millilitre
mm	millimetre
m/z	mass-to-charge ratio
N	Normal; gram molecular weight of solute divided by hydrogen equivalent of solute, per litre of solution
mg	milligram 10^{-3} g
μ g	microgram 10^{-6} g
ng	nanogram 10^{-9} g
pg	picogram 10^{-12} g
fg	femtogram 10^{-15} g
ppm	parts per million (mg/L, mg/kg)
ppb	parts per billion (μ g/L, μ g/kg)
ppt	parts per trillion (ng/L, ng/kg)
ppq	parts per quadrillion (pg/L, pg/kg)
v/v	volume per unit volume
w/v	weight per unit volume
DCM	Dichloromethane (Methylene Chloride)
PFK	Perfluorokerosene
HIRES	High Resolution
GC	Gas Chromatography

MS Mass Spectrometry
HRMS High Resolution Mass Spectrometry

Acceptance Criteria

Values used by the laboratory in order to determine that a process is in control.

Accuracy It is the degree of agreement of a measured value with the true or expected value of the quantity of concern.

Analyte A Nitrosodimethylamine and/or 1,4-Dioxane parameter tested by a method.

Blind Sample It is a sample submitted for analysis whose composition is known to the submitter but unknown to the analyst. A blind sample is used to test the proficiency of a measurement process.

Calibration Standard (CAL)

Consist of a set of solutions containing known amounts of native & carbon-13-labelled NDMA and/or 1,4-Dioxane. These solutions are used to establish the relationship between the parameter's concentration & MS detector response over the expected range of sample concentration.

Calibration Verification Material

Consists of a calibration standard solution of intermediate level concentration (e.g. CS3), used to assess whether the initial calibration is still valid.

Certified Reference Material

It is a stable, homogenous, and well characterized reference material, one or more of whose property values are certified by repetitive analysis by several operators & different methodologies in one or more qualified laboratories of known precision & accuracy. This material is used to assess the accuracy of a measurement process.

CAS# Chemical Abstracts Compound Registry Number.

Control Sample

It is a reference material of known composition that is analyzed concurrently with test samples to evaluate the accuracy and/or precision of a measurement process.

EDL Estimated detection limit or detection limit.

Glassware Proof Rinse

It is the composite final solvent rinse of each piece of glassware intended for use in processing a batch of samples. Proof rinse samples are analyzed before sample processing begins.

Instrument Detection Limit

It is the smallest concentration/amount of analyte, in a solution containing only the analyte(s) of interest, which produces an instrumental response that satisfies all analyte detection & identification criteria.

IS

Internal Standard, a deuterated or ^{13}C -labelled analyte that is added to a sample extract prior to instrument analysis.

Isomer

A member of a group of compounds that differ from each other only in terms of locations of a specified number of common substituent atoms, or groups of atoms, on the parent compound.

Method Blank Laboratory control sample using reagents, purified water, soil or relevant matrix known to be free of contaminants.

Method Detection Limit (MDL)

It is the smallest test sample concentration/amount of analyte that produces an instrumental response that satisfies all analyte detection & identification criteria when the sample is processed & analyzed according to the requirements of a specific test method. Reported MDL values reflect the composite effect of sample-related variables as well as method-related variables.

MSDS

Material Safety Data Sheet

NIOSH

National Institute of Occupational Safety & Health

Precision

It is the degree of agreement between the data generated from repetitive measurements under specified conditions. It is generally reported as the standard deviation (SD) or relative standard deviation (RSD).

%D

Percent Difference.

Quality Assurance (QA)

It is a system of activities whose purpose is to provide the producer or user of a product with the assurance that the product meets a defined standard of quality. The system consists of two separate but related activities, quality control & quality assessment.

Quality Control (QC)

It is the overall system of activities whose purpose is to control the quality of a product so that it meets the needs of users.

Recovery Standards

They are selected compounds that are added to sample extracts immediately before instrumental analysis so that surrogate (internal standard) recoveries can be calculated.

RPD (%) Relative Percent Difference.

Relative Retention Factor (RRF)

It is the quotient of a target analyte response factor (instrument response per unit weight) divided by the response factor (RF) for its corresponding labelled surrogate. An RRF value remains constant over the range of concentration for which instrument response is linear.

RSD Relative Standard Deviation.

SDS Soxhlet/Dean-Stark extractor, an extraction device applied to the extraction of solid & semi-solid materials.

Spiked blank Laboratory control sample that has been fortified with native analytes of interest.

Stock Solution A solution containing an analyte that is prepared using a reference material traceable to EPA, the National Institute of Science & Technology (NIST), or a source that will attest to the purity & authenticity of the reference material.

Surrogate A compound whose composition and chemical properties are nearly identical to those of target analytes, but which is distinguishable from target analytes by some means of detection (i.e. MS). These include deuterated or ¹³C-labelled analogues of the target analytes, which are added to the sample prior to extraction or clean-up steps.

Window Defining Mixture

It is a solution containing the earliest & latest eluting congeners within each homologous group of target analytes on a specified GC column.

SAMPLE DATA

000008

MW-13

Lab Name Maxxam Analytics Inc.

Matrix (soil/water): water

Sample wt/vol: 980 (g/mL) mL

Level (low/med) low

% Moisture Not applicable Decanted (Y/N): N

Concentrated Extract Volume 1000 (uL)

Injection Volume 2 (uL)

Acid Wash Cleanup (Y/N): N pH Not analyzed

Lab Sample ID: A203024-709615

Project Number: 1400

Project Name: JPL

Lab File ID: kr0444

Date Received: February 1, 2002

Date Extracted: February 4, 2002

Lab Batch: 20020204

Date Analyzed: February 4, 2002

Calib. Ref.: 20020131

Time Analyzed: 18:10

Dilution Factor: 1

CAS No.	Compound	Conc. (ug/L)	Qualifier	EDL (ug/L)	RL (ug/L)
62-75-9	NDMA	0.00200	U	0.000390	0.00200
	Surrogate	Recovery (%)	Acceptance Criteria (%)		
76523-40-5	D6-NDMA	24	10-85		

000009

MW-16

Lab Name Maxxam Analytics Inc.

Matrix (soil/water): water

Sample wt/vol: 1000 (g/mL) mL

Level (low/med) low

% Moisture Not applicable Decanted (Y/N): N

Concentrated Extract Volume 1000 (uL)

Injection Volume 2 (uL)

Acid Wash Cleanup (Y/N): N pH Not analyzed

Lab Sample ID: A203024-709618

Project Number: 1400

Project Name: JPL

Lab File ID: kr0447

Date Received: February 1, 2002

Date Extracted: February 4, 2002

Lab Batch: 20020204

Date Analyzed: February 4, 2002

Calib. Ref.: 20020131

Time Analyzed: 18:47

Dilution Factor: 1

CAS No.	Compound	Conc. (ug/L)	Qualifier	EDL (ug/L)	RL (ug/L)
62-75-9	NDMA	0.00200	U	0.000390	0.00200
	Surrogate	Recovery (%)	Acceptance Criteria (%)		
76523-40-5	D6-NDMA	22	10-85		

000015

MW-16D

Lab Name Maxxam Analytics Inc.

Matrix (soil/water): water

Sample wt/vol: 1010 (g/mL) mL

Level (low/med) low

% Moisture Not applicable Decanted (Y/N): N

Concentrated Extract Volume 1000 (uL)

Injection Volume 2 (uL)

Acid Wash Cleanup (Y/N): N pH Not analyzed

Lab Sample ID: A203024-709619

Project Number: 1400

Project Name: JPL

Lab File ID: kr0448

Date Received: February 1, 2002

Date Extracted: February 4, 2002

Lab Batch: 20020204

Date Analyzed: February 4, 2002

Calib. Ref.: 20020131

Time Analyzed: 18:58

Dilution Factor: 1

CAS No.	Compound	Conc. (ug/L)	Qualifier	EDL (ug/L)	RL (ug/L)
62-75-9	NDMA	0.00200	U	0.000390	0.00200
	Surrogate	Recovery (%)	Acceptance Criteria (%)		
76523-40-5	D6-NDMA	20	10-85		

000017

MW-13-MATRIX SPIKE DUPLICATE

Lab Name	Maxxam Analytics Inc.			Lab Sample ID:	A203024-709616
Matrix (soil/water):	water			Project Number	1400
Sample wt/vol:	990	(g/mL)	mL	Project Name:	JPL
Level (low/med)	low			Lab File ID:	kr0445
% Moisture	Not applicable	Decanted (Y/N):	N	Date Received:	February 1, 2002
Concentrated Extract Volume	1000	(uL)		Date Extracted:	February 4, 2002
Injection Volume	2	(uL)		Lab Batch	20020204
Acid Wash Cleanup (Y/N):	N	pH	NA	Date Analyzed:	February 4, 2002
				Calib. Ref.	20020131
				Time Analyzed:	18:22
				Dilution Factor:	1

CAS No.	Compound	MSD Extract Conc (ng/mL)	Extract Conc. (ng/mL)	Spike Level (ng/mL)	Recovery (%)	%RPD MS/MSD	Acceptance Criteria (%)
62-75-9	NDMA	3.92	0.00200	5.00	78	1	25
Internal Standard		Recovery (%)	Acceptance Criteria (%)				
76523-40-5	D6-NDMA	24	10-85				

000013

MW-13-MATRIX SPIKE

Lab Name	Maxxam Analytics Inc.			Lab Sample ID:	A203024-709617
Matrix (soil/water):	water			Project Number	1400
Sample wt/vol:	1000	(g/mL)	mL	Project Name:	JPL
Level (low/med)	low			Lab File ID:	kr0446
% Moisture	Not applicable	Decanted (Y/N):	N	Date Received:	February 1, 2002
Concentrated Extract Volume	1000	(uL)		Date Extracted:	February 4, 2002
Injection Volume	2	(uL)		Lab Batch	20020204
Acid Wash Cleanup (Y/N):	N	pH	NA	Date Analyzed:	February 4, 2002
				Calib. Ref.	20020131
				Time Analyzed:	18:34
				Dilution Factor:	1

CAS No.	Compound	MS Extract Conc. (ng/mL)	Sample Conc. (ng/mL)	Spike Level (ng/mL)	Recovery (%)	Acceptance Criteria (%)
62-75-9	NDMA	3.97	0.00200	5.00	79	27-165
Surrogate		Recovery (%)	Acceptance Criteria (%)			
76523-40-5	D6-NDMA	29	10-85			

000011

Applied P & Ch Laboratory

13760 Magnolia Ave. Chino CA 91710

Tel: (909) 590-1828 Fax: (909) 590-1498

APCL Analytical Report

Submitted to:

SOTA Environmental

Attention: Yu Zeng

16835 W. Bernardo Dr, Ste. 212

San Diego CA 92127

Tel: (858) 485-8100 Fax: (858) 485-0812

Service ID #: 801-021150

Collected by:

Collected on: 01/10/02

Received: 01/10/02

Extracted: N/A

Tested: N/A

Reported: 03/07/02

Sample Description: Water

Project Description: 00HW019 JPL

Analysis of Water Samples

Component Analyzed	Method	Unit	PQL	Analysis Result
				MW-17-3
				02-01150-1

NITROSAMINES BY HRMS ^(a)

PQL: Practical Quantitation Limit. MDL: Method Detection Limit. CRDL: Contract Required Detection Limit

N.D.: Not Detected or less than the practical quantitation limit.

"-": Analysis is not required.


J: Reported between PQL and MDL.

† All results are reported on dry basis for soil samples.

Listed Dilution Factors (DF) are relative to the method default DF. All unlisted DFs are 1.0

^(a) Subcontracted to Maxxam Analytics Inc. See attached.

Respectfully submitted,


Dominic Lau

Laboratory Director

Applied P & Ch Laboratory

Level D Data Package Deliverables

General Information

Project: 010HW019 JPL

APCL Service ID: 02-1150



Applied P & Ch Laboratory
13760 Magnolia Ave. Chino, CA 91710
Telephone (909)590-1828
Fax (909)590-1498



Applied R & Uh Laboratory

APCL
13760 Magnolia Ave. Chino CA 91710
Tel: (909) 590-1828 Fax: (909) 590-1498

Chain of Custody

Please Print in pen Page 2 of 2

Client: SOTA Envir. Tech., Inc.

Contact: Mike Sayre

Tel #: 858-485-8100 Fax #: 858-485-0812

Address: 16835 W. Bernardo Dr. #212

City: San Diego

State: CA

Zip code: 92127

Bill to: SOTA

Project Name/Code SPL

Job # 0040017 P.O. #

Project Address PASADENA, CALIFORNIA APCL Quotation #

Due Date: Regular Rush: days hours Sampled by: MES/SAT

Field Sample ID No.	Sample Description	Date Time Collected	Sample Matrix	Preservation	# of Containers	Analysis Items	White - With report Yellow - Lab copy Pink - Originator
MW-17-1	MW-17-1	1/10/02 1555	WATER	---	1	As (200.9) Cd (200.7) Cr (200.8) Mn (200.9) Ni (200.9) Pb (200.9) Se (200.9) V (200.9) Zn (200.9) Zr (200.9)	
ER-2	EQUIP. RINSE	1425		HCl	3	As (200.9) Cd (200.7) Cr (200.8) Mn (200.9) Ni (200.9) Pb (200.9) Se (200.9) V (200.9) Zn (200.9) Zr (200.9)	
MW-17-3	MW-17-3	1147		---	2	As (200.9) Cd (200.7) Cr (200.8) Mn (200.9) Ni (200.9) Pb (200.9) Se (200.9) V (200.9) Zn (200.9) Zr (200.9)	

1150

QC Requirement: ☐ Regular; ☐ QA/QC Report; ☐ WIP; ☐ Raw Data; ☐ Extended Raw Data ☐ CLP; ☐ ACE ☐ AFCEE ☐ NEBSA (E, C or D); ☒ Other (Please specify)
Sample Disposal: ☐ Return ☐ Disposal by APCL ☐ Hold for days after receiving date. If not specified, samples will be discarded 45 days after samples are received.
Sample Conditions: ☐ Intact; ☐ Broken. Cooler Seal: ☐ Intact; ☐ Broken; ☐ None. Tag # Temperature: ☐ Room ☐ Cold (°C).

Relinquished by *MS* Date/Time 1/10/02 1656 Received by *Relayd Hansen* Date/Time 01-10-02 1656
Relinquished by Date/Time / Received by Date/Time /

APCL USE ONLY Service # Note:

Clients understand that all terms described in the proposals, quotations for this project, and/or the general terms provided in the current APCL price schedules will be followed. APCL reserves the right to terminate its service or withhold delivery of any reports, if in APCL's sole discretion the terms of the project have been broken.

Applied P & Ch Laboratory

13760 Magnolia Ave., Chino CA 91710

Tel: (909) 590-1828 Fax: (909) 590-1828

1150

Sample Receiving Checklist

APCL ServiceID: _____ Client Name/Project: Sota Environmental

1. Sample Arrival

Date/Time Received 1/10/02 656 Date/Time Opened 1/11/02 083 By (name): Paul
Custody Transfer: ☐ Client ☐ Golden State ☐ UPS ☐ US Mail ☐ FedEx ☒ APCL Empl: Richard

2. Chain-of-Custody (CoC)

☒ With Samples? ☐ Faxed? ☒ Client has Copy? ☐ Signed, dated? By: _____
☒ Project ID? ☒ Analyses Clear? ☐ Hold Samples? #on Hold _____ # Received _____
☒ CoC/Docs Zip-Locked under lid? ☐ Compos. #: _____ #Samples OK? _____
☐ Discrepancies? ☐ Client notified? ☐ Response (attach docs): _____

3. Shipping Container/Cooler

☒ Cooler Used? # of _____ Cooled by: ☒ Ice ☐ Blue Ice ☐ Dry Ice ☐ None
Temp °C 2-5
(Cooler temperature measured from temp blank if present, otherwise measured from the cooler).
Cooler Custody Seal? ☐ Absent ☐ Intact ☐ Tampered?

4. Sample Preservation

☐ pH <2 ☐ pH >12
If Not, pH = _____ Preserved by: ☐ Client ☐ APCL ☐ Third Party _____

5. Holding-time Requirements

☐ pH 24hr ☐ BACT 6/24hr ☐ Cr^{VI} 24hr ☐ NO₃⁻ 48hr ☐ BOD 48hr
☐ Cl₂ ASAP ☐ Turbidity 48hr ☐ DO ASAP ☐ Fe(II) ASAP
☐ HT Expired? ☐ Client notified?

6. Sample Container Condition

☒ Intact? ☐ Broken? ☐ Documented? Number: _____
Type: ☒ plastic ☒ glass ☐ Tube: brass/SS ☐ Tedlar Bag
☒ Quantity OK? ☐ Leaking? ☐ Anomaly?
☒ Caps tight? ☐ Air Bubbles? ☐ Anomaly?
Labels: ☐ Unique ID? ☐ Date/Time ☐ Preserved?

7. Turn Around Time

☒ RUSH TAT: 5 day ☐ Std (7-10 days) ☐ Not Marked

8. Sample Matrix

☐ Drinking H₂O ☒ Other Liq ☐ Soil ☐ Wipe ☐ Polymer ☐ Air ☐ Other: _____
☐ Ground H₂O ☐ Sludge ☐ Filter ☐ Oil/Petro ☐ Paint ☐ W. Water ☐ Extract ☐ Unknown

9. Pre-Login Check List Completed & OK?

☒ ALL OK? (if not, attach docs) ☐ Client Contact? (Name: _____) Date/Time: _____
Received/Checked by: Paul Date: 10 Jan 2002 Time: 8:42 a.m.

*HT: Samples must be analyzed for results to reflect total concentrations. Results generated outside required of holding times are considered minimal values and may be used to define waste as hazardous but not as non-hazardous.

DocumentFile: [uncal.textfiles]smprcl.tex.

059300

Applied P & Ch Laboratory

13760 Magnolia Ave. Chino CA 91710

Tel: (909) 590-1828 Fax: (909) 590-1498

Sample Login: Check List

02-01150 (1288_ 298) (4858100_ 298)

01/14/02

Part 1: General Information

<input type="checkbox"/> Company Information	Name:	<i>SOTA Environmental</i>
	Address:	<i>16835 W. Bernardo Dr, Ste. 212 ,San Diego ,CA 92127</i>
<input type="checkbox"/> Project Information	Project Description:	<i>JPL</i>
	Project #:	<i>00HW019</i>
<input type="checkbox"/> Billing Information	P.O. #:	
	Bill Address:	<i>16835 W. Bernardo Dr, Ste. 212 ,San Diego ,CA 92127</i>
	Lab Project ID:	<i>2002.0002</i>
	Client Database #:	<i>01</i>
<input type="checkbox"/> Receiving Information	Who Received Sample?	<i>Paul</i>
	Receiving Date/Time:	<i>01/10/02 1656</i>
	COC No.	
<input type="checkbox"/> Shipping Information	Shipping Company	<i>APCL pick up</i>
	Packing Information:	<i>Cooler/Ice Chester</i>
	Cooler Temperature:	<i>2.5 °C</i>
<input type="checkbox"/> Container Information	Container Provider:	<i>Client</i>
<input type="checkbox"/> Sampling Information	Sampling Person:	
	Sampling Company:	<i>Client</i>
<input type="checkbox"/> Turn-Around-Time Option:		<i>Rush 5 working day(s)</i>
<input type="checkbox"/> QC Option:		<i>NEESA C</i>
<input type="checkbox"/> Disposal Option:		<i>Not specify</i>

059301

Part 2: Sample Information

Seq. #	Sample ID (on COC)	Sample Sub-ID	APCL Sample ID	Matrix	Cont- tainer	Preser- vative	Vol, ml Am. g	# of Replica	Condition G, L, B	Collected mmddyy	Hold ?	Composite Group	TAT Days
1	MW-17-3 / NDMA	02-01150-1	W	G			1000	2	G	011002	N	0	7 <input type="checkbox"/>

Part 3: Analysis Information

Test Items:



Customized-13, Sub-contract

Seq. #	Client's Sample ID (as given on COC)	Sample Sub-ID	APCL Sample ID	Matrix	CUSTOM	
1	MW-17-3	NDMA	02-01150-1	W	X ~	<input type="checkbox"/>

Login By En-Yu Paul Kou

Check By 

Nitrosamines by HRMS



**SUBMISSION CASE NARRATIVE
NDMA**

MAXXAM L.I.M.S. No. A201163

PROJECT: Applied P & Ch Laboratory NDMA Analysis

I. Receipt

Sample was received at Maxxam on January 15, 2002.
Sample was received in good condition.

II. Holding Times

- A. Sample preparation: all holding times were met.
- B. Sample analysis: all holding times were met.

III. Method

The method followed was Maxxam's in-house method for NDMA analysis, Entitled "EXTRACTION & ANALYSIS OF NITROSAMINES AND NDMA BY HRMS" SOP # TO.1021.04.

IV. Preparation

Sample preparation proceeded normally. Sample was extracted on January 16, 2002.

V. Analysis


Analysis proceeded normally. Sample was analyzed on January 21, 2002.

- A. Calibration: All criteria were met.
- B. Mass Resolution: All criteria met.

000001

- C. Method Blank: All acceptance criteria were met for method blank.
- D. Laboratory Control Spike: A LCS and LCSDUP were analyzed and they had a RPD of 10%.
- E. Matrix spike/Matrix spike duplicate: MS and MSD were not analyzed with this sample set.
- F. Surrogate Standards: All samples and QC samples met surrogate Standard criteria
- G. Samples: Sample analysis proceeded normally.
- H. Glass blank: The surrogate recovery on the glass blank exceeds the control limits for water samples. This procedure does not expose the surrogate to any matrix other than the solvent.

I certify that this data package is in compliance with the terms and conditions agreed to by the client and Maxxam Analytics Inc., both technically and for completeness, except for any conditions noted above. Release of the data contained in this hardcopy data package has been authorized by the HRMS Facility Co-ordinator, as verified by the following signature.

 FOR

Patrick Pond
Senior Technical Specialist

SUMMARY OF SAMPLES SUBMITTED-NDMA

(YYYY/MM/DD)

SOTA SAMPLE NO.	MAXXAM L.I.M.S. ID	DATE SAMPLED	DATE RECEIVED	DATE EXTRACTED	DATE ANALYZED	ASSOCIATED QC LABEL
MW-17-3	A201163-698941	2002/01/10	2002/01/15	2002/01/16	2002/01/21	20020116

000003

Glossary of Definitions

NDMA	N-Nitrosodimethylamine
OPR	Ongoing Performance & Recovery Standard (Matrix spike)
PAR	Performance & Recovery Standard (Spiking Mixture)
IPR	Initial Performance & Recovery Standard (Matrix spike)
K-D	Kuderna-Danish concentrator; a device used to concentrate the analytes in a solvent
LIMS	Laboratory Information Management System
MISA	Municipal Industrial Strategy for Abatement
EPA	see USEPA
USEPA	United States Environmental Protection Agency
CEPA	Canadian Environmental Protection Agency
amp	ampere
cm	centimetre
g	gram
h	hour
ID	internal diameter
OD	outside diameter
In.	inch
L	litre
M	Molecular ion
min	minute
mL	millilitre
mm	millimetre
m/z	mass-to-charge ratio
N	Normal; gram molecular weight of solute divided by hydrogen equivalent of solute, per litre of solution
mg	milligram 10^{-3} g
μ g	microgram 10^{-6} g
ng	nanogram 10^{-9} g
pg	picogram 10^{-12} g
fg	femtogram 10^{-15} g
ppm	parts per million (mg/L, mg/kg)
ppb	parts per billion (μ g/L, μ g/kg)
ppt	parts per trillion (ng/L, ng/kg)
ppq	parts per quadrillion (pg/L, pg/kg)
v/v	volume per unit volume
w/v	weight per unit volume
DCM	Dichloromethane (Methylene Chloride)
PFK	Perfluorokerosene
HIRES	High Resolution
GC	Gas Chromatography

MS Mass Spectrometry
HRMS High Resolution Mass Spectrometry

Acceptance Criteria

Values used by the laboratory in order to determine that a process is in control.

Accuracy It is the degree of agreement of a measured value with the true or expected value of the quantity of concern.

Analyte A Nitrosodimethylamine and/or 1,4-Dioxane parameter tested by a method.

Blind Sample It is a sample submitted for analysis whose composition is known to the submitter but unknown to the analyst. A blind sample is used to test the proficiency of a measurement process.

Calibration Standard (CAL)

Consist of a set of solutions containing known amounts of native & carbon-13-labelled NDMA and/or 1,4-Dioxane. These solutions are used to establish the relationship between the parameter's concentration & MS detector response over the expected range of sample concentration.

Calibration Verification Material

Consists of a calibration standard solution of intermediate level concentration (e.g. CS3), used to assess whether the initial calibration is still valid.

Certified Reference Material

It is a stable, homogenous, and well characterized reference material, one or more of whose property values are certified by repetitive analysis by several operators & different methodologies in one or more qualified laboratories of known precision & accuracy. This material is used to assess the accuracy of a measurement process.

CAS# Chemical Abstracts Compound Registry Number.

Control Sample

It is a reference material of known composition that is analyzed concurrently with test samples to evaluate the accuracy and/or precision of a measurement process.

EDL Estimated detection limit or detection limit.

Glassware Proof Rinse

It is the composite final solvent rinse of each piece of glassware intended for use in processing a batch of samples. Proof rinse samples are analyzed before sample processing begins.

Instrument Detection Limit

It is the smallest concentration/amount of analyte, in a solution containing only the analyte(s) of interest, which produces an instrumental response that satisfies all analyte detection & identification criteria.

IS

Internal Standard, a deuterated or ^{13}C -labelled analyte that is added to a sample extract prior to instrument analysis.

Isomer

A member of a group of compounds that differ from each other only in terms of locations of a specified number of common substituent atoms, or groups of atoms, on the parent compound.

Method Blank

Laboratory control sample using reagents, purified water, soil or relevant matrix known to be free of contaminants.

Method Detection Limit (MDL)

It is the smallest test sample concentration/amount of analyte that produces an instrumental response that satisfies all analyte detection & identification criteria when the sample is processed & analyzed according to the requirements of a specific test method. Reported MDL values reflect the composite effect of sample-related variables as well as method-related variables.

MSDS

Material Safety Data Sheet

NIOSH

National Institute of Occupational Safety & Health

Precision

It is the degree of agreement between the data generated from repetitive measurements under specified conditions. It is generally reported as the standard deviation (SD) or relative standard deviation (RSD).

%D

Percent Difference.

Quality Assurance (QA)

It is a system of activities whose purpose is to provide the producer or user of a product with the assurance that the product meets a defined standard of quality. The system consists of two separate but related activities, quality control & quality assessment.

Quality Control (QC)

It is the overall system of activities whose purpose is to control the quality of a product so that it meets the needs of users.

Recovery Standards

They are selected compounds that are added to sample extracts immediately before instrumental analysis so that surrogate (internal standard) recoveries can be calculated.

RPD (%) Relative Percent Difference.

Relative Retention Factor (RRF)

It is the quotient of a target analyte response factor (instrument response per unit weight) divided by the response factor (RF) for its corresponding labelled surrogate. An RRF value remains constant over the range of concentration for which instrument response is linear.

RSD Relative Standard Deviation.

SDS Soxhlet/Dean-Stark extractor, an extraction device applied to the extraction of solid & semi-solid materials.

Spiked blank Laboratory control sample that has been fortified with native analytes of interest.

Stock Solution A solution containing an analyte that is prepared using a reference material traceable to EPA, the National Institute of Science & Technology (NIST), or a source that will attest to the purity & authenticity of the reference material.

Surrogate A compound whose composition and chemical properties are nearly identical to those of target analytes, but which is distinguishable from target analytes by some means of detection (i.e. MS). These include deuterated or ¹³C-labelled analogues of the target analytes, which are added to the sample prior to extraction or clean-up steps.

Window Defining Mixture

It is a solution containing the earliest & latest eluting congeners within each homologous group of target analytes on a specified GC column.

SAMPLE DATA

000008

MW-17-3

Lab Name Maxxam Analytics Inc.

Matrix (soil/water): water

Sample wt/vol: 960 (g/mL) mL

Level (low/med) low

% Moisture Not applicable Decanted (Y/N): N

Concentrated Extract Volume 1000 (uL)

Injection Volume 2 (uL)

Acid Wash Cleanup (Y/N): N pH Not analyzed

Lab Sample ID: A201163-698941

Project Number: 1150

Project Name: JPL

Lab File ID: kr0136

Date Received: January 15, 2002

Date Extracted: January 16, 2002

Lab Batch: 20020116

Date Analyzed: January 21, 2002

Calib. Ref.: 20020118

Time Analyzed: 13:38

Dilution Factor: 1

CAS No.	Compound	Conc. (ug/L)	Qualifier	EDL (ug/L)	RL (ug/L)
62-75-9	NDMA	0.00200	U	0.000390	0.00200
	Surrogate	Recovery (%)	Acceptance Criteria (%)		
76523-40-5	D6-NDMA	39	10-85		

000009

Applied P & Ch Laboratory

13760 Magnolia Ave. Chino CA 91710

Tel: (909) 590-1828 Fax: (909) 590-1498

APCL Analytical Report

Submitted to:

SOTA Environmental

Attention: Yu Zeng

16835 W. Bernardo Dr, Ste. 212

San Diego CA 92127

Tel: (858) 485-8100 Fax: (858) 485-0812

Service ID #: 801-021338

Collected by:

Collected on: 01/24/02

Received: 01/24/02

Extracted: N/A

Tested: N/A

Reported: 02/27/02

Sample Description: Water

Project Description: 00HW019 JPL

Analysis of Water Samples

Component Analyzed	Method	Unit	PQL	Analysis Result
				MW-24-1
				02-01338-1

NITROSAMINES BY HRMS ^(a)

PQL: Practical Quantitation Limit. MDL: Method Detection Limit. CRDL: Contract Required Detection Limit

N.D.: Not Detected or less than the practical quantitation limit.


"-": Analysis is not required.

J: Reported between PQL and MDL.

Listed Dilution Factors (DF) are relative to the method default DF. All unlisted DFs are 1.0

^(a) Subcontracted to Maxxam Analytics Inc. See attached.

Respectfully submitted,



Dominic Lau

Laboratory Director

Applied P & Ch Laboratory

**SUBMISSION CASE NARRATIVE
NDMA**

MAXXAM L.I.M.S. No. A202819

PROJECT: Applied P & Ch Laboratory NDMA Analysis

I. Receipt

Sample was received at Maxxam on January 30, 2002.
Sample was received in good condition.

II. Holding Times

- A. Sample preparation: all holding times were met.
- B. Sample analysis: all holding times were met.

III. Method

The method followed was Maxxam's in-house method for NDMA analysis, Entitled "EXTRACTION & ANALYSIS OF NITROSAMINES AND NDMA BY HRMS" SOP # TO.1021.04.

IV. Preparation

Sample preparation proceeded normally. Sample was extracted on January 31, 2002.

V. Analysis


Analysis proceeded normally. Sample was analyzed on February 4, 2002.

- A. Calibration: All criteria were met.
- B. Mass Resolution: All criteria met.

000001

- C. Method Blank: All acceptance criteria were met for method blank.
- D. Laboratory Control Spike: A LCS and LCSDUP were analyzed and they had a RPD of 7%.
- E. Matrix spike/Matrix spike duplicate: MS and MSD were not analyzed with these samples.
- F. Surrogate Standards: All samples and QC samples met surrogate Standard criteria
- G. Samples: Sample analysis proceeded normally.
- H. Glass blank: All acceptance criteria were met for the glass blank.

I certify that this data package is in compliance with the terms and conditions agreed to by the client and Maxxam Analytics Inc., both technically and for completeness, except for any conditions noted above. Release of the data contained in this hardcopy data package has been authorized by the HRMS Facility Co-ordinator, as verified by the following signature.



AnnMarie Wright
Laboratory Operations Manager

SUMMARY OF SAMPLES SUBMITTED-NDMA									
(YYYY/MM/DD)									
SOTA SAMPLE NO.	MAXXAM	DATE	DATE	DATE	DATE	DATE	DATE	DATE	ASSOCIATED
	L.I.M.S. ID	SAMPLED	RECEIVED	EXTRACTED	ANALYZED	QC LABEL			
MW-24-1	A202819-708292	2002/01/24	2002/01/30	2002/01/31	2002/02/04				20020131

000003

Glossary of Definitions

NDMA	N-Nitrosodimethylamine
OPR	Ongoing Performance & Recovery Standard (Matrix spike)
PAR	Performance & Recovery Standard (Spiking Mixture)
IPR	Initial Performance & Recovery Standard (Matrix spike)
K-D	Kuderna-Danish concentrator; a device used to concentrate the analytes in a solvent
LIMS	Laboratory Information Management System
MISA	Municipal Industrial Strategy for Abatement
EPA	see USEPA
USEPA	United States Environmental Protection Agency
CEPA	Canadian Environmental Protection Agency
amp	ampere
cm	centimetre
g	gram
h	hour
ID	internal diameter
OD	outside diameter
In.	inch
L	litre
M	Molecular ion
min	minute
mL	millilitre
mm	millimetre
m/z	mass-to-charge ratio
N	Normal; gram molecular weight of solute divided by hydrogen equivalent of solute, per litre of solution
mg	milligram 10^{-3} g
µg	microgram 10^{-6} g
ng	nanogram 10^{-9} g
pg	picogram 10^{-12} g
fg	femtogram 10^{-15} g
ppm	parts per million (mg/L, mg/kg)
ppb	parts per billion (µg/L, µg/kg)
ppt	parts per trillion (ng/L, ng/kg)
ppq	parts per quadrillion (pg/L, pg/kg)
v/v	volume per unit volume
w/v	weight per unit volume
DCM	Dichloromethane (Methylene Chloride)
PFK	Perfluorokerosene
HIRES	High Resolution
GC	Gas Chromatography

MS Mass Spectrometry
HRMS High Resolution Mass Spectrometry

Acceptance Criteria

Values used by the laboratory in order to determine that a process is in control.

Accuracy It is the degree of agreement of a measured value with the true or expected value of the quantity of concern.

Analyte A Nitrosodimethylamine and/or 1,4-Dioxane parameter tested by a method.

Blind Sample It is a sample submitted for analysis whose composition is known to the submitter but unknown to the analyst. A blind sample is used to test the proficiency of a measurement process.

Calibration Standard (CAL)

Consist of a set of solutions containing known amounts of native & carbon-13-labelled NDMA and/or 1,4-Dioxane. These solutions are used to establish the relationship between the parameter's concentration & MS detector response over the expected range of sample concentration.

Calibration Verification Material

Consists of a calibration standard solution of intermediate level concentration (e.g. CS3), used to assess whether the initial calibration is still valid.

Certified Reference Material

It is a stable, homogenous, and well characterized reference material, one or more of whose property values are certified by repetitive analysis by several operators & different methodologies in one or more qualified laboratories of known precision & accuracy. This material is used to assess the accuracy of a measurement process.

CAS# Chemical Abstracts Compound Registry Number.

Control Sample

It is a reference material of known composition that is analyzed concurrently with test samples to evaluate the accuracy and/or precision of a measurement process.

EDL Estimated detection limit or detection limit.

Glassware Proof Rinse

It is the composite final solvent rinse of each piece of glassware intended for use in processing a batch of samples. Proof rinse samples are analyzed before sample processing begins.

Instrument Detection Limit

It is the smallest concentration/amount of analyte, in a solution containing only the analyte(s) of interest, which produces an instrumental response that satisfies all analyte detection & identification criteria.

IS

Internal Standard, a deuterated or ^{13}C -labelled analyte that is added to a sample extract prior to instrument analysis.

Isomer

A member of a group of compounds that differ from each other only in terms of locations of a specified number of common substituent atoms, or groups of atoms, on the parent compound.

Method Blank Laboratory control sample using reagents, purified water, soil or relevant matrix known to be free of contaminants.

Method Detection Limit (MDL)

It is the smallest test sample concentration/amount of analyte that produces an instrumental response that satisfies all analyte detection & identification criteria when the sample is processed & analyzed according to the requirements of a specific test method. Reported MDL values reflect the composite effect of sample-related variables as well as method-related variables.

MSDS

Material Safety Data Sheet

NIOSH

National Institute of Occupational Safety & Health

Precision

It is the degree of agreement between the data generated from repetitive measurements under specified conditions. It is generally reported as the standard deviation (SD) or relative standard deviation (RSD).

%D

Percent Difference.

Quality Assurance (QA)

It is a system of activities whose purpose is to provide the producer or user of a product with the assurance that the product meets a defined standard of quality. The system consists of two separate but related activities, quality control & quality assessment.

Quality Control (QC)

It is the overall system of activities whose purpose is to control the quality of a product so that it meets the needs of users.

Recovery Standards

They are selected compounds that are added to sample extracts immediately before instrumental analysis so that surrogate (internal standard) recoveries can be calculated.

RPD (%) Relative Percent Difference.

Relative Retention Factor (RRF)

It is the quotient of a target analyte response factor (instrument response per unit weight) divided by the response factor (RF) for its corresponding labelled surrogate. An RRF value remains constant over the range of concentration for which instrument response is linear.

RSD Relative Standard Deviation.

SDS Soxhlet/Dean-Stark extractor, an extraction device applied to the extraction of solid & semi-solid materials.

Spiked blank Laboratory control sample that has been fortified with native analytes of interest.

Stock Solution A solution containing an analyte that is prepared using a reference material traceable to EPA, the National Institute of Science & Technology (NIST), or a source that will attest to the purity & authenticity of the reference material.

Surrogate A compound whose composition and chemical properties are nearly identical to those of target analytes, but which is distinguishable from target analytes by some means of detection (i.e. MS). These include deuterated or ¹³C-labelled analogues of the target analytes, which are added to the sample prior to extraction or clean-up steps.

Window Defining Mixture

It is a solution containing the earliest & latest eluting congeners within each homologous group of target analytes on a specified GC column.

SAMPLE DATA

000008

MW-24-1

Lab Name Maxxam Analytics Inc.

Matrix (soil/water): water

Sample wt/vol: 970 (g/mL) mL

Level (low/med) low

% Moisture Not applicable Decanted (Y/N): N

Concentrated Extract Volume 1000 (uL)

Injection Volume 2 (uL)

Acid Wash Cleanup (Y/N): N pH Not analyzed

Lab Sample ID: A202819-708292

Project Number: 1338

Project Name: JPL

Lab File ID: kr0422

Date Received: January 30, 2002

Date Extracted: January 31, 2002

Lab Batch: 20020131

Date Analyzed: February 4, 2002

Calib. Ref.: 20020131

Time Analyzed: 12:50

Dilution Factor: 1

CAS No.	Compound	Conc. (ug/L)	Qualifier	EDL (ug/L)	RL (ug/L)
62-75-9	NDMA	0.000910	J	0.000390	0.00200
	Surrogate	Recovery (%)	Acceptance Criteria (%)		
76523-40-5	D6-NDMA	24	10-85		

000009

Applied P & Ch Laboratory

13760 Magnolia Ave. Chino CA 91710

Tel: (909) 590-1828 Fax: (909) 590-1498

APCL Analytical Report

Submitted to:

SOTA Environmental

Attention: Yu Zeng

16835 W. Bernardo Dr, Ste. 212

San Diego CA 92127

Tel: (858) 485-8100 Fax: (858) 485-0812

Service ID #: 801-021729

Collected by:

Collected on: 02/22/02

Sample Description: Water

Project Description: 00HW019 JPL

Received: 02/22/02

Extracted: N/A

Tested: N/A

Reported: 03/22/02

Analysis of Water Samples

Component Analyzed	Method	Unit	PQL	Analysis Result
				MW-7 02-01729-1

NITROSAMINES BY HRMS ^(a)

PQL: Practical Quantitation Limit. MDL: Method Detection Limit. CRDL: Contract Required Detection Limit

N.D.: Not Detected or less than the practical quantitation limit.

"-": Analysis is not required.


J: Reported between PQL and MDL.

† All results are reported on dry basis for soil samples.

Listed Dilution Factors (DF) are relative to the method default DF. All unlisted DFs are 1.0

^(a) Subcontracted to Maxxam Analytics Inc. See attached.

Respectfully submitted,



Dominic Lau

Laboratory Director

Applied P & Ch Laboratory

